(109×210mm Size)

Code No. 299-54501

**Wako**

**Limulus PS Single Test Wako**

**20 Test**

### Instruction


#### [Preface]

Endotoxin is lipopolysaccharide derived from cell membranes of Gram-negative bacteria and is the representative pyrogen. The administration of endotoxin contaminated blood, saline solution and injectants into the human body may induce serious effects such as fever and acute shock. Therefore, it is imperative for therapeutics to be tested for endotoxin contamination. Since reports about coagulation of Limulus amebocyte lysate (LAL) induced by endotoxin derived from bacteria by F. B. Bang in 1965<sup>1)</sup> and J. Levin & F. B. Bang in 1964<sup>2)</sup>, detection of endotoxin using LAL has been widely used as a rapid and sensitive method. The method measuring endotoxin with LAL is already filed by US Pharmacopoeia in 1980, and by Japanese Paemacopeia in 1988<sup>3)4)</sup>.

It has become obvious that LAL is reactive to some limited substances besides bacterial glycolipid. Kakinuma reported that LAL reacts with  $\beta$ -1,3-glucan<sup>5)</sup> and Iwanaga verified a system activating LAL coagulation with  $\beta$ -1,3-glucan apart from that with endotoxin<sup>6)7)</sup>. Moreover, there are reports that some active substance to LAL is eluted from cellulose-based membranes<sup>8)</sup>. All of the examples described above indicate that a special concern must be taken in specificity of LAL to endotoxin. LAL ES reagent used in the test is prepared as endotoxin specific LAL reagent by inactivating  $\beta$ 1,3-glucan induced cascade by adding an excess of a  $\beta$ 1,3-glucan derivative.

Pyrosep used in the test is an affinity resin specific to endotoxin, which is composed of water-insoluble support and histidine as a ligand conjugated through a spacer. This resin, developed by Tanabe Pharmaceuticals, Ltd, is useful for removal of endotoxin from macromolecule solutions and complicated solutions<sup>9)10)</sup>.

#### ⟨Schematic figure of Pyrosep⟩

Condition of absorb : pH 3 ~ 8 ,  $\mu \leq 0.1$   
Capacity of absorb : 0.31mg (*E. coli* UKT-B, LPS)/g absorbant  
Dissociation constant (Kd) :  $7.3 \times 10^{-3}$  M (*E. coli* UKT-B, LPS)

This test is composed of LAL ES reagent and Pyrosep suspension described above. Affinity concentration of endotoxin from sample is done with Pyrosep resin columnchromatography in capillary as the first step, followed by measurement of endotoxin by time resolved turbidometric assay using LAL ES reagent, which allows measurement of endotoxin in much smaller amount than conventional methods allow, even in fatty samples.

#### ⟨Note⟩

This test is not permitted as a replacement for pyrogen test using rabbits.


### [Characteristics]

1. The test measures very small amounts of endotoxin in complicated samples such as kidney artificial dialysis solution without pronounced interference with endogenous materials in the samples, due to affinity concentration of endotoxin by Pyrosep before LAL reaction.
2. The test can be applied to fatty materials such as fat-soluble vitamins if the material is soluble in ethanol because Pyrosep works as affinity resin to endotoxin even in ethanol.
3. The test is cost-conscious because LAL reagent is efficiently used, depending on the numbers of samples to be assayed.

### [Principle]


The coagulation mechanism of LAL is summarized in the following figure. Endotoxin triggers the first step of serine proteinase cascade in LAL. Eventually, coagulogen, at the very end of the cascade, is hydrolyzed into coagulin which forms a soluble gel.  $\beta$ -1,3-glucan also initiates gelation of LAL through another cascade in dose-dependent manner, and high concentrations of glucan inhibits LAL activation without pronounced interference to endotoxin-initiating cascade. LAL ES reagent provided in the test is prepared as a LAL reagent specifically reactive to endotoxin by adding excess  $\beta$ -1,3-glucan derivative<sup>1)</sup>.

Time resolved turbidimetric assay of endotoxin with a Toxinometer stoichiometrically determines the concentration of endotoxin in samples based on measurement of gelation time (Tg) between turbidity changes from initiation of the reaction to the threshold of the predetermined transmittance level.

Coagulation mechanism of LAL

In the first step of the test, endotoxin in the sample is specifically concentrated by passing the sample through Pyrosep resin packed in glass capillary column. Absorbed endotoxin on the resin is directly conducted to LAL ES reagent and the coagulation level specifically induced by endotoxin is measured with a Toxinometer.

### [Content]

1. LAL ES reagent, *Limulus polyphemus* amebocyte lysate lyophilized (containing Tris-HCl buffer and  $\beta$ -1,3-glucan derivative) ..... 20 vials for 0.3 mL
2. Pyrosep resin suspension (containing in phosphate buffer) ..... 20 vials 0.77 mL
3. LAL reconstitution solution ..... 2 vials ( 7 mL)
4. Wash solution ..... 4 vials (11 mL)
5. Sample diluent ..... 4 vials (10 mL)
6. Glass capillary column (having a filter at the connected end with silicon tube, endotoxin free) ..... 20 pieces

### (Necessary materials and instruments to be prepared)

- 1) a syringe with 10-20 mL capacity (without needle) ..... 1 piece
- 2) a Stopcock with a 3 ways ..... 1 piece
- 3) a syringe fitting for 1.5 mmID tube ..... 1 piece
- 4) a size 10 binder clip (about 51 mm binding site)  
as indicated in the below figure 4 ..... 1 piece
- 5) Toxinometer ..... 1 instrument


### [Assay procedure]

#### I. Preparation of assay sample


Using sample diluent, adjust the pH of the sample to about 5. Generally, pH in most of the samples and endotoxin solutions for the standard curve can be adjusted, simply by adding sample diluent with 1/10 volume of the sample. If buffer action of the sample is relatively strong, concentrated acid solution can be used for pH adjustment. The final sample volume should not be increased significantly for the subsequent experiment.

#### II. Procedure

1. Open Pyrosep resin suspension vial by removing aluminum cap. Suspend Pyrosep resin in the vial by several shaking and inverting.
2. Pick up a glass capillary attached with silicon tube at their connection site. Do not touch glass capillary part after removal of aluminum foil, as shown in figure 1.
3. Place the free end of the glass capillary into the suspension, and connect the free end of silicon tune to a syringe via a 3 way stopcock and a syringe fitting prepared as shown in figure 2.
4. Gently agitate the suspension and draw all of the suspension into the capillary column by drawing the pestle of the syringe. Then switch the stopcock to close the channel connecting to the capillary and to open another channel to drain the solution in the syringe, as shown in figure 3. For this step, Pyrosep resin is packed in the capillary column by stacking the resin to the filter integrated inside the capillary.

(Fig 1)

(Fig 2)

- Carefully transfer and place the end of the glass capillary into the sample solution in the sample tube without touching the glass part. Extract all of the solution after turning the stopcock to open the capillary channel as indicated in figure 3.

**⟨Notes⟩**

- A binder clip, as indicated in figure 4, is useful to fix the plunger in place, as illustrated in figure 4.
- It takes about 20 minutes to transfer 5mL of sample to the column.
- Discard the solution in the syringe by switching the stopcock to drain and pushing the plunger.

**⟨Note⟩**

The sample applied resin column should be washed with 1-2 mL of wash solution provided in the test if the sample is feared to contain some factors activating or inhibiting gelation of LAL.

In the same manner as described in 5 and 6.

- Reconstitute LAL ES reagent with 0.3 mL of reconstitution solution provided.

**⟨Note⟩**


Preferably, LAL ES reagent should be reconstituted just before use. It should be used within 1 hour after reconstitution.

- Transfer the capillary end to the LAL reagent, and draw it slowly by drawing the plunger, while observing the solution rising inside the capillary.
- Turn the stopcock off the channel connecting to the capillary just when the LAL reagent reaches the position at about 6 cm beneath the capillary-silicon tube connection site, and intrude the air through another channel into the syringe cavity. Then, return the stopcock to connect the syringe and the capillary, pushing back the resin and the LAL with the air into the LAL ES vial.

**⟨Notes⟩**

- On extruding the resin from the capillary column, the free end of the capillary should be placed in the position of the vial inner wall, 1 cm above the surface of the residual LAL reagent in order to avoid creating bubbles in the LAL reagent, as shown in figure 5.


(Fig 3)

- ▶ When suck of Pyrosep resin suspension, sample and reconstitute LAL ES, → Close ②
- ▶ When it drain the solution in the syringe or intrude the air into the syringe, → Close ①


(Fig 4)

How to use of a file binder

(Fig 5)

How to push back the LAL

- b) Abrupt extrusion of the resin and LAL reagent from the capillary may cause the filter integrated inside the capillary to slip out.
- c) Immediate extrusion of the LAL reagent and the resin is required after turning the stopcock to close the capillary. Incomplete switching from the LAL reagent extraction to the air intrusion may cause continued rising of the LAL reagent in the capillary, causing it to reach the syringe fitting which may be contaminated with endotoxin.
- d) A small quantity of the resin remaining in the capillary will not cause problems in results of the assay. LAL extraction and extrusion should be repeated if a large quantity of the resin is retained in the capillary.
10. Check the LAL-resin suspension whether or not the suspension contains bubbles, and eliminate bubbles of more than 2 mm in diameter by tapping the vial wall with finger if necessary.
11. After vortexing the vial for about 5 seconds, apply it to Toxinometer to determine gelation time<sup>12)13)14)15)</sup>.

#### ⟨Note⟩

Generally, the threshold value is determined at the level of transmittance of light at 94.9%. (Wait time 5.0 min.)


12. Calibrate the concentration of endotoxin in the sample using the standard curve plotted with value obtained from measurement of endotoxin solutions reference in the same manners as described above.

#### ⟨Notes⟩

The standard curve of endotoxin is plotted as X bar: log(Endotoxin concentration);Y bar: log[log (gelation time)], to determine the concentration with linear regression. Refer to the leaflet, "Toxinometer Standard Assay Procedure" for details.

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