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Wako

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
1. BAMBANKER® (LYMPHOTEC Inc.)
2. Medium with serum (含血清, Company A)
3. Serum-free Medium (无血清, Company A)

human γ T cells

Storage Period	10 months *1
Bambanker	65%
Serum Medium	37%
Serum-free Medium	35%


Cell number/vial 1.0×10^6

human B cell line

Storage Period	9 months *1
Bambanker	74%
Serum Medium	54%
Serum-free Medium	35%


Cell number/vial 1.0×10^6

PC12 (rat-derived adrenal pheochromocytoma)

Storage Period	11 months *1
Bambanker	95%
Serum Medium	59%
Serum-free Medium	20%


Cell number/vial 1.0×10^6

P3U1 (mouse myeloma cell line)

Storage Period	1 year *1
Bambanker	95%
Serum Medium	95%
Serum-free Medium	70%


Cell number/vial 2.0×10^6

OKT4 (mouse hybridoma)

Storage Period	1 year *1
Bambanker	100%
Serum Medium	100%
Serum-free Medium	92%


Cell number/vial 1.3×10^6

human gastric epithelial cells

Storage Period	10 months *1
Bambanker	100%
Serum Medium	62%
Serum-free Medium	56%


Cell number/vial 1.0×10^6

monkey B cell line

Storage Period	10 months *1
Bambanker	56%
Serum Medium	40%
Serum-free Medium	18%

Cell number/vial 1.0×10^6

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- Zhang C., Seo J., Nakamura T. (2018) Cellular Approaches in Investigating Argonaute2-Dependent RNA Silencing. In: Okamura K., Nakanishi K. (eds) Argonaute Proteins. Methods in Molecular Biology, vol 1680. Humana Press, New York, NY.
- Sharma, A., Mücke, M., & Seidman, C. E. (2018). Human induced pluripotent stem cell production and expansion from blood using a non-integrating viral reprogramming vector. Current Protocols in Molecular Biology, 122, e58. doi: 10.1002/cpm.58.
- Souta Motoike, Mikihiro Kajiya, Nao Komatsu, et al. Cryopreserved clumps of mesenchymal stem cell/extracellular matrix complexes retain osteogenic capacity and induce bone regeneration. Stem Cell Res Ther. 2018; 9: 73. Published online 2018 Mar 21. doi: 10.1186/s13287-018-0826-0.
- Srijaya Thekkeparambil Chandrabose, Sandhya Sriram, et al. Amenable epigenetic traits of dental pulp stem cells underlie high capability of xeno-free episomal reprogramming. Stem Cell Research & Therapy 2018 9:68.
- Evans, Michael A. et al. "Macrophage-Mediated Delivery of Light Activated Nitric Oxide Prodrugs with Spatial, Temporal and Concentration Control." Chemical Science (2018): n. pag. Web. doi:10.1039/C8SC00015H.
- Burrige P W, Diecke S, Matsa E, et al. Modeling cardiovascular diseases with patient-specific human pluripotent stem cell-derived cardiomyocytes [J]. Patient-Specific Induced Pluripotent Stem Cell Models: Generation and Characterization, 2016: 119-130
- Durruthy-Durruthy J, Sebastiano V, Wossidlo M, et al. The primate-specific noncoding RNA HPAT5 regulates pluripotency during human preimplantation development and nuclear reprogramming[J]. Nature genetics, 2016, 48(1): 44-52.
- Kojima S I, Kaku M, Kawata T, et al. Cranial suture-like gap and bone regeneration after transplantation of cryopreserved MSCs by use of a programmed freezer with magnetic field in rats[J]. Cryobiology, 2015, 70(3): 262-268.

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FUJIFILM Wako Pure Chemical Corporation