CiMSTM Instruction for use



Product Description

CiMS[™] is a medium for culture of human mesenchymal stem cells (hMSCs). CiMS[™] is optimized for the multiple passage expansion of hMSCs (isolated from Bone Marrow, Adipose Tissue and Dental Pulp) while still allowing for differentiation into the several desired lineages. Xeno-free or Animal-free supplement can be selected according to the intended use. In addition, there is no need to pre-coat culture vessels with any type of attachment matrix before using CiMS[™].

Product	Catalog Number (NIPRO/CSTI)	Components	Volume	Storage
CiMS [™] -BM	87-070 / A2G00P05C	Basal media (Animal-free)	500 mL	2-8 ℃ ; Protect from Light
*CiMS [™] -sXF	87-071 / A2G10P1CC	Xeno-free supplement	10 mL	-20 ℃ ; Protect from Light
*CiMS TM -sAF	87-072 / A2G20P1CC	Animal-free supplement	10 mL	-20 ℃ ; Protect from Light
Related Product	Catalog Number (NIPRO/CSTI)	Components	Volume	Storage
rTE	87-974 / 1210	Recombinant Trypsin/EDTA solution	100 mL	2-8 $^\circ \!$
sTI	87-975 / 1220	Synthetic Trypsin inhibitor	100 mL	2-8 $^\circ \!$
PBS(-)	87-949, 87-972 / 1102P05, 1102P10		500, 1000 mL	1-30 ℃

* Please select either Xeno-free or Animal-free supplement. Xeno-free : It contains human derived component. Any other animal derived component free. Animal-free : animal derived component free

Storage instructions

 CiMS[™] instructions: upon arrival, store CiMS[™]-BM(basal medium) protected from light at 2°C to 8°C and CiMS[™]-sXF or CiMS[™]-sAF at ≤ -20°C in a freezer. Once thawed, supplement should be stored 2-8°C and added to basal medium within 2 weeks. After supplement is added to basal medium, use immediately.

Preparation of Culture Media

- Decontaminate the external surfaces of the CiMS[™]-BM and CiMS[™]-sXF or −sAF bottle with 70% v/v ethanol or isopropanol.
- Add 1/50 volume of CiMS[™]-sXF or CiMS[™]-sAF to CiMS[™]-BM (CiMS[™] complete medium).
- * Recommend to make necessary volume of the medium just before use.

Initiation of Culture Process

- Rapidly thaw a frozen cryotube of hMSC in a 37℃ water bath until a small amount of ice remains.
- 2. Wipe cryotube with 70% v/v ethanol or isopropanol before opening. Pipet the entire contents of cryotube into a conical tube.
- Carefully add 10 mL of pre-warmed (37°C) CiMS[™] complete medium.
- **4.** Centrifuge the tubes at about100 x g for 5 minutes at room temperature and discard the supernatant.
- 5. Resuspend the cell pellet in pre-warm (37°C) CiMS[™] complete medium and add the cell suspension to an appropriate culture vessel (cell culture vessels for adherent cells) at a density of 0.5 – 1.0 x 10⁴ cells / cm².

6. Incubate at 37℃, 5% CO₂, humidified incubator.

Maintenance

1. Change CiMS[™] complete medium every 3 days.

Subculturing

The following instructions are for a 25 cm² flask. Adjust all volumes accordingly for other size vessels.

- **1.** Subculture the cells when they are about 85% confluent.
- **2.** Remove the culture medium from 25 cm² flask.
- **3.** Cover the cells with 0.5 mL rTE (Catalog #87-974).
- Place the culture vessel into a 37℃ humidified incubator for 3-5 minutes. Periodically examine the cell layer microscopically and check for cell detachment.
- **5.** Allow the trypsinization to continue until approximately 90% of the cells are rounded up.
- **6.** At this point, tap the flask gently to release the majority of cells from the culture surface. If only a few cells detach, you may not have let them trypsinize long enough.
- 7. After cell released, neutralize the trypsin in the vessel with 0.5 mL of s-TI (Catalog #87-975) at room temperature.
- **8.** Quickly transfer the detached cells to sterile 15 mL centrifuge tube.

- **9.** Rinse the flask with a final 10 mL of PBS(-)(Catalog #87-949, 87-972) to collect residual cells, and add this rinse to the centrifuge tube.
- **10.** Examine the harvested flask under microscope to make sure the harvest was successful by looking at the number of cells left behind. This should be less than 5%.
- **11.** Centrifuge the tubes at about 100 x g for 5 minutes at room temperature and discard the supernatant.
- 12. Resuspend the pellet in a minimum volume of temperature equilibrated CiMS[™] complete medium by gently pipetting up and down. Count total number of viable cells.
- **13.** Add the calculated volume of cell suspension to each flask and gently rock to disperse the cell suspension over the growth surface.
- **14.** Incubate at 37°C, 5% CO₂, humidified incubator.

Products are for research use only and not intended for human or animal diagnostic or therapeutic uses unless otherwise stated.



CELL SCIENCE & TECHNOLOGY INSTITUTE INC.

CD17372



北京 Tel: 010 64136388 上海 Tel: 021 62884751 广州 Tel: 020 87326381 香港 Tel: 852 27999019 info@boppard.cn www.boppard.cn





产品数据卤