A hand is shown holding a nautilus shell, which is a natural example of a logarithmic spiral. The shell is held between the thumb and index finger, with the rest of the hand visible in the lower left. The background is a dark, textured grey. The text 'MicroRNA Research Products' is overlaid on the shell in a bold, orange font.

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## References

1. Bartel DP: MicroRNA:genomics, biogenesis, mechanism, and function. *Cell* 2004, 116:281-297.
2. Cortez MA, Calin GA : MicroRNA identification in plasma and serum: a new tool to diagnose and monitor diseases. *Expert Opin Biol Ther* 2009, 9:703-711
3. Hammond SM, Bernstein E, Beach D, Hannon GJ: An RNA-directed nuclease mediates post-transcriptional gene silencing in *Drosophila* cells. *Nature* 2000, 404:293-296.
4. Hammond SM, Boettcher S, Caudy AA, Kobayashi R, Hannon GJ: Argonaute2, a link between genetic and biochemical analyses of RNAi. *Science* 2001, 293:1146-1150.
5. Gregory RI, Chendrimada TP, Cooch N, Shiekhattar R: Human RISC couples microRNA biogenesis and posttranscriptional gene silencing. *Cell* 2005, 123:631-640.
6. Sasaki T, Shiohama A, Minoshima S, Shimizu N: Identification of eight members of the Argonaute family in the human genome. *Genomics* 2003, 82:323-330.
7. Pillai RS, Artus CG., Filipowicz W: Tethering of human Ago proteins to mRNA mimics the miRNA-mediated repression of protein synthesis. *RNA* 2004, 10:1518-1525.
8. Meister G, Landthaler M, Patkaniowska A, Dorsett Y, Teng G, Tuschl T: Human Argonaute2 mediates RNA cleavage targeted by miRNAs and siRNAs. *Mol Cell* 2004, 15:185-197.
9. Liu J, Carmell MA, Rivas FV, Marsden CG, Thomson JM, Song JJ, Hammond SM, Joshua-Tor L, Hannon GJ: Argonaute2 is the catalytic engine of mammalian RNAi. *Science* 2004, 305:1437-1441.
10. Song JJ, Smith SK, Hannon GJ, Joshua-Tor L: Crystal structure of Argonaute and its implications for RISC slicer activity. *Science* 2004, 305:1434-1437.
11. Meister G, Tuschl T: Mechanisms of gene silencing by double-stranded RNA. *Nature* 2004, 431:343-349.
12. Hutvagner G, Zamore PD: A microRNA in a multiple-turnover RNAi enzyme complex. *Science* 2002, 297:2056-2060.
13. Yekta S, Shih IH, Bartel DP: MicroRNA-directed cleavage of HOXB8 mRNA. *Science* 2004, 304:594-596.
14. Beitzinger M, Peters L, Zhu JY, Kremmer E, Meister G: Identification of human microRNA targets from isolated Argonaute protein complexes. *RNA Biol* 2007, 4:76-84.
15. Karginov FV, Conaco C, Xuan Z, Schmidt BH, Parker JS, Mandel G, Hannon GJ: A biochemical approach to identifying microRNA targets. *Proc Natl Acad Sci USA* 2007, 104:19291-19296.
16. Easow G, Teleman AA, Cohen SM: Isolation of microRNA targets by miRNP immunopurification. *RNA* 2007, 13:1198-1204.
17. Hendrickson DG, Hogan DJ, Herschlag D, Ferrell JE, Brown PO: Systematic identification of mRNAs recruited to Argonaute2 by specific microRNAs and corresponding changes in transcript abundance. *PloS ONE* 2008, 3:e2126.
18. Landthaler M, Gaidatzis D, Rothballer A, Chen PY, Soll SJ, Dinic L, Ojo T, Hafner M, Zavolan M, Tuschl T: Molecular characterization of human Argonaute-containing ribonucleoprotein complexes and their bound target mRNAs. *RNA* 2008, 14:2580-2596.
19. Chi SW, Zang JB, Mele A, Darnell RB: Argonaute HITS-CLIP decodes microRNA-mRNA interaction maps. *Nature* 2009, 460:479-486.
20. Wang WX, Wilfred BR, Hu Y, Stromberg AJ, Nelson PT: Anti-Argonaute RIP-Chip shows that miRNA transfections alter global patterns of mRNA recruitment to microribonucleoprotein complexes. *RNA* 2009 Dec.30.
21. Hayashida Y, Nishibu T, Inoue K, Kurokawa T: A useful approach to total analysis of RISC associated RNA . *BMC Res Notes* 2009, 2:169.

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1. Anti-Human Ago2 Magnetic Beads Solution 600μL×1瓶
2. Cell Lysis Solution 20mL×1瓶
3. Washing Solution for IP 40mL×1瓶
4. Elution Solution I for IP 500μL×1瓶
5. Elution Solution II for IP 500μL×1瓶
6. Binding Solution for Spin Column 2mL×1瓶
7. Binding Enhancer for Spin Column 100μL×1瓶
8. Washing Solution I for Spin Column 3mL×1瓶
9. Washing Solution II for Spin Column 4mL×1瓶
10. Elution Solution for Spin Column 1mL×1瓶
11. Spin Column/Collection Tube 10瓶

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### hAgo2(4G8) MAb

Cell. 2014 Jan 16;156:146-157.

Inefficient SRP interaction with a nascent chain triggers a mRNA quality control pathway.

### microRNA Isolation Kit human / mouse Ago1, and human Ago2

Nucleic Acids Research , 2014,1-13.

Novel functional small RNAs are selectively loaded onto mammalian Ago1.

### microRNA Cloning Kit Wako, ssDNA Ligase, thermostable

Mol Cell. 2011 Nov 4;44(3):424-36.

MCPIP1 Ribonuclease Antagonizes Dicer and Terminates MicroRNA Biogenesis through Precursor MicroRNA Degradation.

Suzuki HI, Arase M, Matsuyama H, Choi YL, Ueno T, Mano H, Sugimoto K, Miyazono K.

### Ago1(2A7), hAgo2(4G8), hAgo3(1C12) MAb

RNA Biol. 2011 Jan 1;8(1).

Deep-sequencing of human Argonaute-associated small RNAs provides insight into microRNA sorting and reveals Argonaute association with RNA fragments of diverse origin.

Burroughs AM, Ando Y, Hoon ML, Tomaru Y, Suzuki H, Hayashizaki Y, Daub CO.,

### hAgo2(4G8) MAb

PNAS 2011 Mar 29.

MicroRNA let-7 establishes expression of  $\beta$ 2-adrenergic receptors and dynamically down-regulates agonist-promoted down-regulation.

Wang WC, Juan AH, Panebra A, Liggett SB.

### mAgo2(2D4) MAb

EMBO J . 2011 March 2; 30(5): 823-834.

Small RNA-mediated regulation of iPS cell generation Zhonghan Li,Chao-Shun Yang, Katsuhiko Nakashima, and Tariq M Rana

### mAgo2(2D4) MAb

Molecular Pharmacology March 22, 2011

RISC bound siRNA is a determinant of RNAi mediated gene silencing in mice.

Wei J, Jones J, Kang J, Card A, Krimm M, Hancock P, Pei Y, Ason B, Payson E, Dubinina N, Cancilla M, Stroh M, Burchard J, Sachs A, Hochman J, Flanagan WM, Kuklin N.

### mAgo2(2D4) MAb

Nature structural & molecular biology 2011, 18, 2, Feb

Genome-wide identification of Ago2 binding sites from mouse embryonic stem cells with and without mature microRNA.

Leung AK, Young AG, Bhutkar A, Zheng GX, Bosson AD, Nielsen CB, Sharp PA.

### Ago1(2A7), hAgo2(4G8), hAgo3(1C12) MAb

### microRNA Isolation Kit, Ago1, Human Ago2 and Human Ago3

Genome Res. August 18, 2010

A comprehensive survey of 39 animal microRNA modification events and a possible role for 39 adenylation in modulating microRNA targeting effectiveness.

Burroughs AM, Ando Y, de Hoon MJ, Tomaru Y, Nishibu T, Ukekawa R, Funakoshi T, Kurokawa T, Suzuki H, Hayashizaki Y, Daub CO.

### **hAgo2(4G8) MAb**

**Nucleic Acids Research** , 2009, 1–13

Expanded RNA-binding activities of mammalian Argonaute 2.

Tan GS, Garchow BG, Liu X, Yeung J, Morris JP 4th, Cuellar TL, McManus MT, Kiriakidou M.

### **hAgo2(4G8) MAb**

**Nature structural & molecular biology** 2009 Dec 16;12:1259

Distinct passenger strand and mRNA cleavage activities of human Argonaute proteins.

Wang B, Li S, Qi HH, Chowdhury D, Shi Y, Novina CD

### **microRNA Isolation Kit, Human Ago2**

**BIOLOGY OF REPRODUCTION** 2009 81, 717–729

Human Villous Trophoblasts Express and Secrete Placenta-Specific MicroRNAs into Maternal Circulation via Exosomes.

Luo SS, Ishibashi O, Ishikawa G, Ishikawa T, Katayama A, Mishima T, Takizawa T, Shigihara T, Goto T, Izumi A, Ohkuchi A, Matsubara S, Takeshita T, Takizawa T.

### **hAgo2(4G8) MAb**

**BMC Research Notes** 2009, 2:169

A useful approach to total analysis of RISC-associated RNA.

Hayashida Y, Nishibu T, Inoue K, Kurokawa T.

### **hAgo2(4G8) MAb**

**RNA** 2009, 15:1078–1089

Mammalian GW182 contains multiple Argonaute-binding sites and functions in microRNA-mediated translational repression.

Takimoto K, Wakiyama M, Yokoyama S.

### **mAgo2(2D4) microRNA Isolation Kit, Human Ago2**

**Reproduction** 2008 Dec;136(6):811-22.

MicroRNA (miRNA) cloning analysis reveals sex differences in microRNA expression profiles between adult mouse testis and ovary.

Mishima T, Luo SS, Ishibashi O, Kawahigashi Y, Mizuguchi Y, Ishikawa T, Mori M, Kanda T, Goto T, Takizawa T.

### **hAgo2(4G8) MAb**

**Methods Mol Biol.** 2008;442:29-43.

In vitro RNA cleavage assay for Argonaute-family proteins.

Miyoshi K, Uejima H, Nagami-Okada T, Siomi H, Siomi MC.

### **hAgo2(4G8) MAb**

**Nature** , 2008 Sep 18;455(7211):421-4.

Prolyl 4-hydroxylation regulates Argonaute 2 stability.

Qi HH, Ongusaha PP, Myllyharju J, Cheng D, Pakkanen O, Shi Y, Lee SW, Peng J, Shi Y.

### **hAgo2(4G8) MAb**

**PNAS** 2008 June 10, vol. 105 no. 23 7964–7969

Characterization of endogenous human Argonautes and their microRNA partners in RNA silencing.

Azuma-Mukai A, Oguri H, Mituyama T, Qian ZR, Asai K, Siomi H, Siomi MC.

### **hAgo2(4G8) MAb**

**J Hepatol.** 2009 Mar;50(3):453-60. Epub 2008 Jul 9.

Regulation of the hepatitis C virus genome replication by miR-199a.

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Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit our online catalog to search for other products from Wako; <http://www.e-reagent.com> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

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