

Detection of AP using the Alkaline Phosphatase Staining Kit II

Overview

The following procedure describes staining one well of a 6-well plate for Alkaline Phosphatase (AP) detection.

Product Description	Cat. No.	Format
REPROCELL Stemgent® Alkaline Phosphatase Staining Kit II	00-0055	500 rxn
Components	Size	Storage
Fix Solution	25 mL	4 °C
AP Staining Solution A	10 mL	4 °C
AP Staining Solution B	10 mL	4°C
AP Staining Solution C	10 mL	4 °C

Table 1: REPROCELL Stemgent® Alkaline Phosphatase Staining Kit II contents.

For other plate formats, see the suggested amounts in Table 2.

Additional Materials Required

- 1× PBS
- Tween® 20
- Mounting medium (optional)
- 15 mL conical tubes

Material Preparation

PBST

In a 15 mL conical tube, add 10 mL of 1× PBS. Add 5 μ L of Tween 20 for a final concentration of 0.05%. Mix well and store at room temperature.

• AP Substrate Solution

For one well of a 6-well plate, mix 0.5 mL of Solution A and 0.5 mL of Solution B in a 15 mL conical tube. Incubate at room temperature for 2 minutes. Add 0.5 mL of Solution C.

Note: Prepare only the amount of AP Substrate Solution necessary for each experiment. Quantities can be scaled up or down, as long as a 1:1:1 ratio is preserved. For optimal results, the AP Substrate Solution should be used within 30 minutes after preparation. Discard any remaining solution.

AP Staining of Cells

- 1. Aspirate the culture medium and wash the cells with 2 mL of $1 \times PBST$.
- 2. Add 1 mL of Fix Solution and incubate at room temperature for 2 to 5 minutes.

Note: Do not over fix the cells. Excessive fixation will result in the loss of AP activity.

- 3. Aspirate the Fix Solution and wash the fixed cells with 2 ml of 1× PBST. Do not allow the wells to dry.
- 4. Aspirate the 1× PBST and add 1.5 mL of freshly prepared AP Substrate Solution.
- 5. Incubate the cells in the dark (wrapped with foil or in a dark container) at room temperature for 5 to 15 minutes.

Note: Closely monitor the color change and stop the reaction when the color turns bright to avoid non-specific staining.

- 6. Stop the reaction by aspirating the AP Substrate Solution and washing the wells twice with 2 mL of $1 \times PBS$.
- 7. Cover the cells with 1× PBS or mounting medium to prevent drying.
- 8. AP expression will result in a red or purple stain, while the absence of AP expression will result in no stain.
- 9. Store the plate at 4 °C

Culture vessel	Surface area/well	Fix Solution	1× PBS	AP Staining Solution	Reactions Per kit
24-well plate	2 cm ²	0.5 mL	0.5 mL	0.6 mL	50
12-well plate	3.8 cm ²	1 mL	1 mL	1 mL	24
6-well plate	9.6 cm ²	2 mL	2 mL	1.5 mL	12

Table 2. Suggested amounts per well.

Tween® is a registered trademark of ICI America, Inc

富士胶片和光(广州)贸易有限公司

广州市越秀区先烈中路69号东山广场30楼3002-3003室

北京 Tel: 13611333218 上海 Tel: 021 62884751 广州 Tel: 020 87326381 香港 Tel: 852 27999019

询价: wkgz.info@fujifilm.com 官网: labchem.fujifilm-wako.com.cn





REPROCELL USA Inc

9000 Virginia Manor Road Suite 207 Beltsville, MD 20705 USA

Tel: +1 301 470 3362 Email: info-us@reprocell.com

REPROCELL Europe Ltd

Thomson Pavilion, Todd Campus West of Scotland Science Park Acre Road Glasgow G20 0XA UK

Tel: +44 (0)141 465 3460 Email: <u>info-emea@reprocell.com</u>

REPROCELL India Ltd

3-1-135/1A, CNR Complex Mallapur Hyderabad 500 076 Telangana India

Tel: +44 (0)141 465 3460

Email: Bhargavi.Gurram@reprocell.com

REPROCELL Inc (Japan)

MetLife Shin-yokohama 381, Bldg. 9F 3-8-11, Shin-yokohama Kohokuku, Yokohama Kanagawa 222-0033 Japan

Tel: +81 45 475 3887

Email: info-asia@reprocell.com

www.reprocell.com



PR-STEMGENT-00-0055-D001001