

## Ampdirect 推荐使用酶

### ■ Native and recombinant Taq (天然、重组 Taq 酶)

TaKaRa Taq	(TaKaRa)
Taq DNA Polymerase, Recombinant	(Invitrogen)
AmpliTaq DNA Polymerase	(ABI)
Taq DNA Polymerase	(SIGMA)
Taq DNA Polymerase	(QIAGEN)
Taq DNA Polymerase	(Roche)
Taq DNA Polymerase	(STRATAGENE)
NovaTaq DNA Polymerase	(Novagen)

### ■ Hot Start Enezyme (抗体型热启动酶)

TaKaRa Taq Hot Start Version	(TaKaRa)
Platinum Taq DNA Polymerase	(Invitrogen)
JumpStart Taq DNA Polymerase	(SIGMA)

### ■ Hot Start Enezyme (化学修饰型热启动酶)

NovaTaq Hot Start DNA Polymerase	(Novagen)
----------------------------------	-----------

## Ampdirect 不推荐使用酶

### ■ 除上述推荐以外的化学修饰型热启动酶

AmpliTaq Gold DNA Polymerase	(ABI)
FastStart Taq DNA Polymerase	(Roche)
SureStart Taq DNA Polymerase	(STRATAGENE)
etc.	

### ■ 除了 Taq 之外的 DNA 聚合酶