



3小时, HCP即可现形!

Host Cell Protein ELISAs (ENZ-KIT127, ENZ-KIT128 & ENZ-KIT162)

对CHO, *E. Coli*和HEK293T来源的生物制品中HCP残留污染进行高灵敏检测


生物制药常以细胞培养表达系统来生产和纯化,会产生内源蛋白形式的污染,即宿主细胞蛋白(Host Cell Proteins, HCP)。HCP 杂质可引发免疫反应,因此必须尽量控制生物制品中的污染并在进行监测。Enzo 提供高灵敏的 ELISA 试剂盒,用于检测以 CHO, *E. coli* 和 HEK293T 表达系统生产的生物制品中的 HCP。

特点

- 专为生物制品中 HCP 污染检测而设
- 灵敏度高,检测下限低至 10 ng/mL (CHO), 30 ng/mL (*E. coli*) 或 37 ng/mL (HEK293T)
- 3 小时可以得到高通量检测的结果

产品列表

Sensitive Measurement of HCPs



Typical standard curve for CHO HCP ELISA Kit (ENZ-KIT128)

细胞来源	产品编号	产品名称	规格	介绍
CHO	ENZ-KIT128-0001	CHO Host Cell Protein ELISA Kit	96 wells	灵敏度: 10 ng/mL 检测范围: 3~810 ng/mL
	相关试剂			
	ENZ-PRT122	CHO Host Cell Protein Standards	50 µg	ELISA, WB用 浓度: 1 mg/mL
	ENZ-ABS260	CHO Host Cell Proteins Polyclonal Antibody	100 µg	WB用 免源多抗
	ENZ-ABS261-0100	CHO Host Cell Proteins Polyclonal Antibody (Biotin Conjugate)	100 µg	WB用 免源多抗, 生物素偶联
<i>E. Coli</i>	ENZ-KIT127-0001	<i>E. Coli</i> Host Cell Protein ELISA Kit	96 wells	灵敏度: 30 ng/mL 检测范围: 3~810 ng/mL
	相关试剂			
	ENZ-PRT121	<i>E. Coli</i> Host Cell Protein Standards	50 µg	ELISA, WB用 浓度: 0.5 mg/mL
	ENZ-ABS262	<i>E. Coli</i> Host Cell Proteins Polyclonal Antibody	100 µg	WB用 免源多抗
	ENZ-ABS263-0100	<i>E. Coli</i> Host Cell Proteins Polyclonal Antibody (Biotin Conjugate)	100 µg	WB用 免源多抗, 生物素偶联
HEK293T	ENZ-KIT162-0001	HEK293T Host Cell Protein ELISA Kit 	96 wells	灵敏度: 37 ng/mL 检测范围: 37~27000 ng/mL

富士胶片和光(广州)贸易有限公司

广州市越秀区先烈中路69号东山广场30楼
3002-3003室
询价: info@boppard.cn
目录价查询: bb-china.net

www.boppard.cn

北京 Tel: 010 64136388
上海 Tel: 021 62884751
广州 Tel: 020 87326381
香港 Tel: 852 27999019

官方微信



目录价查询

