

JCSS级别离子标准溶液

概要

JCSS(日本校准服务体系)是用于审查企业是否符合国际标准化机构及国际电力标准会议所规定的校准机关要求事项(ISO/IEC 17025)，并开展校准企业注册的服务体系。

Wako的pH标准溶液、金属标准溶液、离子标准溶液、有机标准溶液(挥发性有机化合物)产品均在JCSS进行注册，溯源性符合值得信赖的国家计量标准，因此能够从根本上保证实验、校准结果的可靠性。


IAJapan(独立行政法人产品评价技术基盘机构认可中心)，已加盟ILAC(国家实验室认可合作组织)及APLAC(亚太实验室认可合作组织)。Wako经IAJapan批准，成为国际MRA对应认可企业。Wako能够发行带有ILAC·MRA及JCSS认可标志的校准证明书。

JCSS证书载的数值通过ILAC/APLAC的MRA认证，能在国际范围内被接受。

*Wako作为校准机关、标准物质的生产者，适用于ISO Guide34:2009及ISO/IEC 17025:2005(JIS Q17025:2005)。(标准物质生产者的认可不在ILAC/MRA范围内)。

特点

能附上带有ILAC MRA及JCSS认可标志的校准证明书



产品一览

离子标准溶液

产品编号	产品名称	规格
019-24011	Anion Mixture Standard Solution 1 阴离子混合标准溶液1	50 mL
010-26361	Anion Mixture Standard Solution 2 阴离子混合标准溶液2	50 mL
019-15461	Ammonium Ion Standard Solution (NH4+ 1000) 铵离子标准溶液 (NH4+ 1000)	50 mL
024-15331	Bromide Ion Standard Solution (Br- 1000) 溴离子标准溶液 (Br- 1000)	50 mL
066-03401	Fluoride Ion Standard Solution (F- 1000) 氟离子标准溶液 (F- 1000)	50 mL

产品编号	产品名称	规格
143-06441	Nitrate Ion Standard Solution (NO ₃ - 1000) 硝酸根离子标准溶液 (NO ₃ - 1000)	50 mL
140-06451	Nitrite Ion Standard Solution (NO ₂ - 1000) 亚硝酸根离子标准溶液 (NO ₂ - 1000)	50 mL
168-17461	Phosphate Ion Standard Solution (PO ₄ (3-) 1000) 磷酸离子标准溶液 (PO ₄ (3-) 1000)	50 mL
192-10821	Sulfate Ion Standard Solution(SO ₄ (2-) 1000) 硫酸根离子标准溶液 (SO ₄ (2-) 1000)	50 mL

有机标准溶液

产品编号	产品名称	规格
221-01971	16 VOC Mixture Standard Solution (Methanol Solution) 16种VOC混合物标准溶液	2 mL×5A
224-01721	23 VOC Mixture Standard Solution (Methanol Solution) 23种VOC混合物标准溶液 (每种物质浓度为1 mg/mL, 甲醇溶液中)	2 mL×5A
066-06561	Formaldehyde Standard Solution (Methanol Solution) 甲醛标准溶液 (甲醇溶液)	2 mL×5A

pH标准溶液

产品编号	产品名称	规格
151-01845	Oxalate pH Standard Solution 草酸盐pH标准溶液	500 mL
166-12141	Phthalate pH Standard Solution	100 mL
168-12145	邻苯二甲酸盐pH标准溶液	500 mL
163-12151	Phosphate pH Standard Equimolar Solution	100 mL
165-12155	中性磷酸盐pH标准溶液	500 mL
166-17445	Phosphate pH Standard Solution 磷酸盐pH标准溶液	500 mL
203-08771	Tetraborate pH Standard Solution	100 mL
205-08775	四硼酸盐pH标准溶液	500 mL
037-16145	Carbonate pH Standard Solution (pH 10.01) 碳酸盐pH标准溶液 (pH 10.01)	500 mL

金属标准溶液

产品编号	产品名称	规格
016-18271	Aluminium Standard Solution (Al 100) 铝标准溶液 (Al 100)	100 mL
016-15471	Aluminium Standard Solution (Al 1000) 铝标准溶液 (Al 1000)	100 mL
013-18281	Antimony Standard Solution (Sb 100) 锑标准溶液 (Sb 100)	100 mL
010-15491	Antimony Standard Solution (Sb 1,000) 锑标准溶液 (Sb 1,000)	100 mL

产品编号	产品名称	规格
013-15501	Arsenic Standard Solution (As 100) 砷标准溶液 (As 100)	100 mL
013-15481	Arsenic Standard Solution (As 1000) 砷标准溶液 (As 1000)	100 mL
027-15321	Barium Standard Solution (Ba 1000) 钡标准溶液 (Ba 1000)	100 mL
023-14201	Bismuth Standard Solution (Bi 100) 铋标准溶液 (Bi 100)	100 mL
021-12661	Bismuth Standard Solution 铋标准溶液 (Bi 1,000)	100 mL
025-16581	Boron Standard Solution(B 1000) 硼标准溶液 (B 1000)	100 mL
030-16211	Cadmium Standard Solution (Cd 100) 镉标准溶液 (Cd 100)	100 mL
036-16171	Cadmium Standard Solution (Cd 1000) 镉标准溶液 (Cd 1000)	100 mL
036-17891	Calcium Standard Solution (Ca 100) 钙标准溶液 (Ca 100)	100 mL
039-16161	Calcium Standard Solution (Ca 1000) 钙标准溶液 (Ca 1000)	100 mL
030-21341	Cesium Standard Solution(Cs 1000) 铯标准溶液 (Cs 1000)	100 mL
037-16221	Chromium Standard Solution (Cr 100) 铬标准溶液 (Cr 100)	100 mL
030-16191	Chromium Standard Solution (Cr 1000) 铬标准溶液 (Cr 1000)	100 mL
039-17901	Cobalt Standard Solution (Co 100) 钴标准溶液 (Co 100)	100 mL
033-16181	Cobalt Standard Solution(Co 1000) 钴标准溶液 (Co 1000)	100 mL
034-16231	Copper Standard Solution (Cu 100) 铜标准溶液 (Cu 100)	100 mL
033-16201	Copper Standard Solution (Cu 1000) 铜标准溶液 (Cu 1000)	100 mL
091-03851	Iron Standard Solution (Fe 100) 铁标准溶液 (Fe 100)	100 mL
094-03841	Iron Standard Solution (Fe 1000) 铁标准溶液 (Fe 1000)	100 mL
070-05781	Gallium Standard Solution(Ga 1000) 镓标准溶液 (Ga 1000)	100 mL
092-05841	Indium Standard Solution(In 1000) 铟标准溶液 (In 1000)	100 mL
127-04301	Lead Standard Solution (Pb 100) 铅标准溶液 (Pb 100)	100 mL
124-04291	Lead Standard Solution (Pb 1000) 铅标准溶液 (Pb 1000)	100 mL
129-05221	Lithium Standard Solution (Li 1000) 锂标准溶液 (Li 1000)	100 mL
136-13601	Magnesium Standard Solution (Mg 100) 镁标准溶液 (Mg 100)	100 mL

产品编号	产品名称	规格
136-12121	Magnesium Standard Solution (Mg 1000) 镁标准溶液 (Mg 1000)	100 mL
139-12111	Manganese Standard Solution (Mn 100) 锰标准溶液 (Mn 100)	100 mL
133-12131	Manganese Standard Solution (Mn 1000) 锰标准溶液 (Mn 1000)	100 mL
135-13671	Mercury Standard Solution (Hg 100) 汞标准溶液 (Hg 100)	100 mL
138-13661	Mercury Standard Solution (Hg 1000) 汞标准溶液 (Hg 1000)	100 mL
130-14961	Molybdenum Standard Solution (Mo 1000) 钼标准溶液 (Mo 1000)	100 mL
144-06471	Nickel Standard Solution (Ni 100) 镍标准溶液 (Ni 100)	100 mL
147-06461	Nickel Standard Solution (Ni 1000) 镍标准溶液 (Ni 1000)	100 mL
162-19941	Potassium Standard Solution (K 100) 钾标准溶液 (K 100)	100 mL
165-17471	Potassium Standard Solution (K 1000) 钾标准溶液 (K 1000)	100 mL
188-01951	Rubidium Standard Solution (Rb 1000) 铷标准溶液 (Rb 1000)	100 mL
192-13861	Selenium Standard Solution (Se 1000) 硒标准溶液 (Se 1000)	100 mL
191-12111	Sodium Standard Solution (Na 100) 钠标准溶液 (Na 100)	100 mL
199-10831	Sodium Standard Solution (Na 1000) 钠标准溶液 (Na 1000)	100 mL
199-13871	Strontium Standard Solution (Sr 1000) 锶标准溶液 (Sr 1000)	100 mL
209-17921	Tellurium Standard Solution (Te 1000) 碲标准溶液 (Te 1000)	100 mL
205-16301	Thallium Standard Solution (Tl 1000) 铊标准溶液 (Tl 1000)	100 mL
202-16311	Tin Standard Solution (Sn 1000) 锡标准溶液 (Sn 1000)	100 mL
221-01851	Vanadium Standard Solution (V 1000) 钒标准溶液 (V 1000)	100 mL
261-01431	Zinc Standard Solution (Zn 100) 锌标准溶液 (Zn 100)	100 mL
264-01421	Zinc Standard Solution (Zn 1000) 锌标准溶液 (Zn 1000)	100 mL

上述试剂仅供实验研究用,不可用作“医药品”、“食品”、“临床诊断”等。

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit our online catalog to search for other products from FUJIFILM Wako; <https://labchem-wako.fujifilm.com> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

富士胶片和光(广州)贸易有限公司

广州市越秀区先烈中路69号东山广场30楼
3002-3003室
询价: info@boppard.cn
目录价查询: bb-china.net

www.boppard.cn

北京 Tel: 010 64136388
上海 Tel: 021 62884751
广州 Tel: 020 87326381
香港 Tel: 852 27999019

化学分析微信



目录价查询

