

日立氨基酸分析仪配套试剂

——缓冲液、显色液、氨基酸标准溶液

日立氨基酸分析仪国内使用的型号主要为L-8500型, L-8500A型, L-8800型, L-8900型等, 其配套试剂分为三部分:

● 缓冲液

● 显色液

● 氨基酸混合标准溶液

缓冲液

日立前期推出的包括L-8500型, L-8500A型, L-8800型氨基酸分析仪配套的缓冲液均分为PH和PF系列, 适用于不同的样本和检测目的。L-8900型可以使用L-8500型的配套试剂。

PH系列: 适用蛋白质水解产物分析

仪器为L-8500型, L-8500A型, L-8800型, 检测样本为蛋白质水解产物。

产品编号	产品名称	包装
634-12471	MCI Buffer L-8500-PH-1	1L×2
604-07463		1L×10
608-07466	MCI Buffer L-8500-PH-2	1L×10
631-12481		1L×2
606-07462	MCI Buffer L-8500-PH-3	1L×10
638-12491		1L×2
631-12501	MCI Buffer L-8500-PH-4	1L×2
602-07464		1L×10
638-12511	MCI Buffer L-8500-PH-RG	1L×2
600-07465		1L×10
608-07461	MCI Buffer L-8500-PH Kit 包含: PH-11L*2, PH-21L*1, PH-31L*1, PH-41L*2, PH-RG1L*1	1kit

PF系列：适用生物液体分析

仪器型号是L-8500型, L-8500A型, L-8800型, 检测样本为生物液体。

产品编号	产品名称	包装
601-07473	MCI Buffer L-8500-PF-1	1L×10
639-12421		1L×2
607-07475	MCI Buffer L-8500-PF-2	1L×10
636-12431		1L×2
633-12441	MCI Buffer L-8500-PF-3	1L×2
605-07476		1L×10
630-12451	MCI Buffer L-8500-PF-4	1L×2
609-07474		1L×10
603-07472	MCI Buffer L-8500-PF-RG	1L×10
637-12461		1L×2
605-07471	MCI Buffer L-8500-PF Kit 包含: PF-11L*2, PF-21L*1, PF-31L*1, PF-41L*2, PF-RG1L*1	1kit

显色液

日立氨基酸分析仪配套的茚三酮显色液 (Ninhydrin Coloring Solution Kit), 适合L-8500, L-8800和L-8900等多种型号日立氨基酸分析仪。

产品编号	产品名称	包装
299-70501	Ninhydrin Coloring Solution Kit for HITACHI 日立氨基酸分析仪配套茚三酮显色液	2L

氨基酸混合标准溶液

Wako生产的氨基酸混合标准溶液, 可用于氨基酸自动分析仪或高效液相色谱(HPLC)进行氨基酸成分分析或测序。

产品编号	产品名称	包装
013-08391	Amino Acids Mixture Standard Solution, Type H 氨基酸混合标准溶液, H型	5mL
015-14461	Amino Acids Mixture Standard Solution, Type AN-2 氨基酸混合标准溶液, AN-2型	5mL
016-08641	Amino Acids Mixture Standard Solution, Type B 氨基酸混合标准溶液, B型	5mL

上述试剂仅供实验研究用, 不可用作“医药品”、“食品”、“临床诊断”等。

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit our online catalog to search for other products from FUJIFILM Wako; <https://labchem-wako.fujifilm.com> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

富士胶片和光(广州)贸易有限公司

广州市越秀区先烈中路69号东山广场30楼
3002-3003室

北京 Tel: 010 64136388

上海 Tel: 021 62884751

广州 Tel: 020 87326381

香港 Tel: 852 27999019

询价: wkqz.info@fujifilm.com

官网: labchem.fujifilm-wako.com.cn

官方微信



目录价查询

