

注射标准用水




我们注射标准用水的质量可满足您所有细胞生长的用途

- 符合生物制药、生命科学和健康行业生产流程的严格标准
 - 细胞培养基和生物试剂的生产
 - 设备的清洁和冲洗
- 质量优于美、欧、日三版药典要求
- 可检测细胞生长的毒性和细胞病变效应
- 配置的预处理系统、气压式蒸馏水器和注射配水回路, 可实现高产
- 适用于各个规格的生物反应器



分步水处理和注射用水生产



可靠的性能

测定项	技术参数
外观	无色透明液体
无菌性	合格
内毒素含量	< 0.25 EU/mL
电导率	5 μ S/cm
可氧化物质含量	符合标准
酸碱度	符合标准
氯化物含量	符合标准
硝酸盐含量	NMT 0.2 ppm
硫酸盐含量	符合标准
TOC 含量	\leq 8000 磅
铵盐含量	NMT 0.2 ppm
钙镁含量	符合标准
颗粒物含量	符合标准

我们的注射标准用水水质符合甚至优于目前美国药典、欧洲药典和日本药典规定的注射用无菌水的测试需求。每个 WFI 批次都已通过美欧日三本药典所规定的包装无菌注射用水标准测试，并配备质量分析证书。

我们的注射标准用水在符合 cGMP 要求的设施内生产，保证最高水准的质量和性能。多种包装和规格可支持从研发到大规模生产的各种需要。

订购信息

标准包装	目录号 #	规格*	备注
WFI 注射标准用水	9309-1 L	1 L	瓶装
	9309-20 LC	20 L	大玻璃瓶装
	9309-20 L	20 L	生物过程容器装；纸箱装运2D端端口袋
	9309-200 L	200 L	生物过程容器装；桶装运2D面端口枕式袋

*可根据要求定制包装尺寸：5 L、10 L、100 L、500 L。

您可以联系 FUJIFILM Irvine Scientific 销售代表或访问我们的官网 www.irvinesci.com 了解关于更多公司与产品的信息。

上述试剂仅供实验研究用，不可用作“医药品”、“食品”、“临床诊断”等。

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit our online catalog to search for other products from FUJIFILM Wako: <https://labchem-wako.fujifilm.com/> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

富士胶片 and 光 (广州) 贸易有限公司

广州市越秀区先烈中路69号东山广场30楼
3002-3003室

北京 Tel: 13611333218

上海 Tel: 021 62884751

广州 Tel: 020 87326381

香港 Tel: 852 27999019

询价: wkgz.info@fujifilm.com

官网: labchem.fujifilm-wako.com.cn

官方微信



目录价查询

