


注射标准用水

我们注射标准用水的质量可满足您所有细胞生长的用途

- 符合生物制药、生命科学和健康行业生产流程的严格标准
 - 细胞培养基和生物试剂的生产
 - 设备的清洁和冲洗
- 质量优于美、欧、日三版药典要求
- 可检测细胞生长的毒性和细胞病变效应
- 配置的预处理系统、气压式蒸馏水器和注射配水回路, 可实现高产
- 适用于各个规格的生物反应器



分步水处理和注射用水生产



可靠的性能

测定项	技术参数
外观	无色透明液体
无菌性	合格
内毒素含量	< 0.25 EU/mL
电导率	5 µS/cm
可氧化物质含量	符合标准
酸碱度	符合标准
氯化物含量	符合标准
硝酸盐含量	NMT 0.2 ppm
硫酸盐含量	符合标准
TOC 含量	≤8000 磅
铵盐含量	NMT 0.2 ppm
钙镁含量	符合标准
颗粒物含量	符合标准

我们的注射标准用水水质符合甚至优于目前美国药典、欧洲药典和日本药典规定的注射用无菌水的测试需求。每个 WFI 批次都已通过美欧日三本药典所规定的包装无菌注射用水标准测试，并配备质量分析证书。

我们的注射标准用水在符合 cGMP 要求的设施内生产，保证最高水准的质量和性能。多种包装和规格可支持从研发到大规模生产的各种需要。

订购信息

标准包装	目录号 #	规格*	备注
WFI注射标准用水	9309-1 L	1 L	瓶装
	9309-20 LC	20 L	大玻璃瓶装
	9309-20 L	20 L	生物过程容器装；纸箱装运2D端端口袋
	9309-200 L	200 L	生物过程容器装；桶装运2D面端口枕式袋

*可根据要求定制包装尺寸：5 L、10 L、100 L、500 L。

上述试剂仅供实验研究用，不可用作“医药品”、“食品”、“临床诊断”等。

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit our online catalog to search for other products from FUJIFILM Wako: <https://labchem-wako.fujifilm.com/> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

富士胶片 and 光 (广州) 贸易有限公司
广州市越秀区先烈中路69号东山广场30楼
3002-3003室

北京 Tel: 13611333218
上海 Tel: 021 62884751
广州 Tel: 020 87326381
香港 Tel: 852 27999019

询价: wkgz.info@fujifilm.com
官网: labchem.fujifilm-wako.com.cn

官方微信



目录价查询



