

体外诊断试剂研发及生产用原料

FUJIFILM Wako提供各类可用于体外诊断试剂研发及生产的酶、表面活性剂、稳定剂、反应辅助因子、显色底物及标准品等高品质试剂, 满足生化、免疫等体外诊断试剂的研发和生产需求。

■ 临床生化检测原料

肝功能

项目	常见方法	名称	产品编号	等级	规格	功能	
ALT/ AST/ LDH	IFCC酶偶联速率法	Tris碱	207-06275	试剂特级	500 g	缓冲体系	
		还原型辅酶 I 二钠	042-16233	生物化学用	500 mg	显色反应	
			046-16231		100 mg		
GLDH	酶速率法	三乙醇胺 (pH 8.0)	142-05615	和光一级	500 mL	缓冲体系	
ALP	IFCC AMP缓冲液法	AMP	014-02586	和光一级	500 mL	缓冲体系	
GGT	IFCC Szasz法	Tris碱	207-06275	试剂特级	500 g	缓冲体系	
TBIL	重氮法	牛磺胆酸钠	199-10032	和光特级	25 g	反应辅助试剂	
DBIL	胆红素氧化酶法	Triton X-100	167-21101	生物化学用	3 L	反应辅助试剂	
			166-11801	化学用	3 L		
			168-11805		500 mL		
TBA	循环酶法	Tris碱	207-06275	试剂特级	500 g	缓冲体系	
ALB	溴甲酚绿 (BCG) 法	溴甲酚绿	020-02682	试剂特级	25 g	显色反应	
PA	免疫比浊法	PEG 6000	169-09125	和光一级	500 g	反应辅助试剂	
			Tween-20	167-11515	化学用	500 g	反应辅助试剂
				166-21115	生物化学用	500 g	

肾功能

项目	常见方法	名称	产品编号	等级	规格	功能
Crea	肌氨酸氧化酶法 (酶法)	4-AAP	011-02275	试剂特级	500 g	显色反应
		TOOS	346-08811	—	1 g	显色反应
Urea	尿素酶-谷氨酸脱氢酶法	脲酶	216-00783	生物化学用	500 mg	酶
			210-00781		1 g	
		Tris碱	207-06275	试剂特级	500 g	缓冲体系
		还原型辅酶 I 二钠	042-16233	生物化学用	500 mg	显色反应
046-16231	100 mg					
UA	尿酸酶-过氧化物酶法	尿酸酶	218-00721	生物化学用	50 U	缓冲体系
		4-AAP	011-02275	试剂特级	500 g	显色反应
		TOOS	346-08811	—	1 g	显色反应

心肌

项目	常见方法	名称	产品编号	等级	规格	功能
CysC	免疫比浊法 (PETIA)	PEG 6000	169-09125	和光一级	500 g	反应辅助试剂
		Tween-20	167-11515	化学用	500 g	反应辅助试剂
			166-21115	生物化学用	500 g	
		Tris碱	207-06275	试剂特级	500 g	缓冲体系
CK	IFCC酶偶联法	咪唑	095-00015	和光特级	500 g	缓冲体系
Hcy	循环酶法	Tris碱	207-06275	试剂特级	500 g	缓冲体系

血脂

项目	常见方法	名称	产品编号	等级	规格	功能
TC / HDL-C / LDL-C	胆固醇氧化酶法 / 直接法 (选择性表面活性剂) / 直接法 (两步法或保护性 表面活性剂)	胆固醇酯酶	033-11223	生物化学用	500U	酶
			037-11221		1000U	
		山葵过氧化物酶	169-10791	生物化学用	100 mg	酶
			165-10793		1 g	
		Triton X-100	167-21101	和光特级	3 L	反应辅助试剂
			166-11801	化学用	3 L	
					168-11805	
		4-AAP	011-02275	试剂特级	500 g	显色反应
TOOS	346-08811	—	1 g	显色反应		
TG	GPO-PAP法	脂蛋白脂肪酶	129-04501	生物化学用	250 mg	酶
		Tween-20	167-11515	化学用	500 g	反应辅助试剂
			166-21115	生物化学用	500 g	
		Triton X-100	167-21101	生物化学用	3 L	反应辅助试剂
			166-11801	化学用	3 L	
					168-11805	
		4-AAP	011-02275	试剂特级	500 g	显色反应
TOOS	346-08811	—	1 g	显色反应		

血糖

项目	常见方法	名称	产品编号	等级	规格	功能
GLU	己糖激酶法	Tris碱	207-06275	试剂特级	500 g	缓冲体系
	葡萄糖氧化酶法	葡萄糖氧化酶	074-02401	生物化学用	20000 U	酶
		4-AAP	011-02275	试剂特级	500 g	显色反应
		TOOS	346-08811	—	1 g	显色反应
HbA1C	酶法 (果糖基氨基酸氧化酶)	Triton X-100	167-21101	生物化学用	3 L	反应辅助试剂
			166-11801	化学用	3 L	
					168-11805	
		4-AAP	011-02275	试剂特级	500 g	显色反应
	TOOS	346-08811	—	1 g	显色反应	
	免疫比浊法	PEG 6000	169-09125	和光一级	500 g	反应辅助试剂
		Tween-20	167-11515	化学用	500 g	反应辅助试剂
			166-21115	生物化学用	500 g	
Tris碱		207-06275	试剂特级	500 g	缓冲体系	

■ 特色产品

产品名称	产品编号	等级	规格	CAS号	用途
脂蛋白脂肪酶	129-04501	生物化学用	250 mg	9004-02-8	酶
海藻糖二水合物	206-18455	和光特级	500 g	6138-23-4	蛋白稳定剂
脱氧胆酸钠 (DOC、SDC)	194-08311	生物化学用	100 g	302-95-4	蛋白增溶剂
酪蛋白酸钠	143-02845	化学用	500 g	9005-46-3	封闭剂、稳定剂
聚乙二醇 6000 (PEG-6000)	169-09125	和光一级	500 g	25322-68-3	增敏剂、促聚剂
	167-22941	分子生物学用	1 kg		分子诊断样本前处理
聚乙烯吡咯烷酮K30	161-03105	—	500 g	9003-39-8	蛋白稳定剂、封闭液、冻干保护剂
	165-17035	和光特级	500 g	9003-39-8	
聚乙烯吡咯烷酮K90	168-03115	—	500 g	9003-39-8	免疫比浊促进剂
	162-17045	和光特级	500 g	9003-39-8	
聚乙二醇 (10) 辛基苯基醚 (Triton X-100)	168-11805	化学用	500 mL	9002-93-1	反应加速剂、反应体系稳定剂、 样本/试剂增溶剂
	167-21101		3 L		
	166-11801	生物化学用	3 L		
咪唑	095-00015	和光特级	500 g	288-32-4	缓冲体系
三羟甲基氨基甲烷 (Tris碱)	207-06275	试剂特级	500 g	77-86-1	缓冲体系
DA-67	046-22341	生物化学用	50 mg	115871-18-6	显色反应

■ 体外诊断试剂原料产品一览

酶

产品名称	产品编号	等级	规格	CAS号	来源
脂蛋白脂肪酶	129-04501	生物化学用	250 mg	9004-02-8	<i>Pseudomonas sp.</i> 铜绿单胞菌属
胆固醇酯酶	033-11223	生物化学用	500 U	9026-00-0	<i>Pseudomonas sp.</i> 铜绿单胞菌属
	037-11221		1000 U		
beta-D-半乳糖苷酶	072-04141	生物化学用	10000 U	9031-11-2	<i>E. Coli</i> 大肠杆菌
胆碱氧化酶	037-14401	生物化学用	500 U	9028-67-5	<i>Alcaligenes sp.</i> 产碱杆菌属(重组)
蛋白酶K, 重组	160-28916	生物化学用	1 g	39450-01-6	<i>Pichia pastoris</i> 毕赤酵母
蛋白酶K溶液, 重组	161-28701	分子生物学用	5 mL	39450-01-6	<i>Tritirachium album</i> 白色念球菌
碱性磷酸酶	018-10693	生物化学用	50 U	9001-78-9	<i>E. Coli</i> 大肠杆菌
抗坏血酸氧化酶	012-10853	生物化学用	200 U	9029-44-1	<i>Cucurbita sp.</i> 南瓜属
亮氨酸脱氢酶, 重组型	120-04391	生物化学用	200 U	9082-71-7	<i>E. Coli</i> 大肠杆菌
尿酸酶	218-00721	生物化学用	50 U	9002-12-4	<i>Candida sp.</i> 念珠菌属
脲酶	216-00783	生物化学用	500 mg	9002-13-5	Jack bean 刀豆/剑豆
	210-00781		1 g		
葡糖淀粉酶(根霉)	073-02851	生物化学用	5000 U	9032-08-0	<i>Rhizopus sp.</i> 根霉属
葡萄糖脱氢酶	071-02411	生物化学用	100 U	9028-53-9	<i>Bacillus sp.</i> 杆菌属
葡萄糖氧化酶	074-02401	生物化学用	20000 U	9001-37-0	<i>Aspergillus niger</i> 黑曲霉
山葵过氧化物酶	162-12883	EIA用	20000 U	9003-99-0	Horseradish 辣根
	165-10793	生物化学用	1 g		
	169-10791		100 mg		
酸性磷酸酶, 土豆来源	010-18193	生物化学用	1000 U	9001-77-8	Potato 土豆
	014-18191		200 U		
	018-18194		50 U		
纤维素酶, 热稳定, 重组溶液	030-19871	生物化学用	1 mL	/	<i>E. Coli</i> 大肠杆菌
脂肪酶	353-46081	—	10 g	9001-62-1	<i>Burkholderia cepacia</i> LP-7 洋葱伯克霍尔德菌

多功能原料

产品名称	产品编号	等级	规格	CAS号
海藻糖二水合物	206-18455	和光特级	500 g	6138-23-4
脱氧胆酸钠 (DOC、SDC)	194-08311	生物化学用	100 g	302-95-4
十二烷基三甲基溴化铵 (DTAB)	359-45525	—	500 g	1119-94-4
油酸钠	194-02635	化学用	500 g	143-19-1
海藻酸钠80-120	196-13325	和光一级	500 g	9005-38-3
硫酸钠	197-03345	试剂特级	500 g	7757-82-6
氯化钠	195-15975	药物测试用	500 g	7647-14-5
柠檬酸钠	203-13605	和光特级	500 g	68-04-2
丁二酸	194-04335	和光特级	500 g	110-15-6
牛磺胆酸钠	199-10032	和光特级	25 g	145-42-6

反应辅助因子

产品名称	产品编号	等级	规格	CAS号
聚乙烯吡咯烷酮K90 PVP K90	162-17045	和光特级	500 g	9003-39-8
	168-03115	—	500 g	
聚乙二醇6,000 (PEG-6,000)	169-09125	和光一级	500 g	25322-68-3
	167-22941	分子生物学用	1 kg	
聚乙二醇4,000 (PEG-4,000)	162-09115	和光一级	500 g	25322-68-3
2-羟乙基-β-环糊精	324-84233	—	100 g	128446-35-5
六水氯化镁, 99.9%	135-09205	—	500 g	7791-18-6
99.9%氯化镁六水合物	139-09203	—	50 g	7791-18-6
六水合氯化镁	131-00162	试剂特级	25 g	7791-18-6
一水乙酸钙	031-00315	试剂特级	500 g	5743-26-0
乙酸镁四水合物	130-00095	和光特级	500 g	16674-78-5
酒石酸钾钠四水合物	192-02991	和光特级	100 g	6381-59-5

表面活性剂

产品名称	产品编号	等级	规格	CAS号
聚氧乙烯(20)失水山梨醇单油酸酯 Tween 80 (吐温80)	164-21775	化学用	500 g	9005-65-6
聚氧乙烯(20)失水山梨醇单棕榈酸酯 Tween 40 (吐温40)	164-11525	化学用	500 g	9005-66-7
聚氧乙烯失水山梨醇月桂酸酯 Tween 20 (吐温20)	166-21115	生物化学用	500 g	9005-64-5
	167-11515	化学用	500 g	
聚乙二醇(10)辛基苯基醚 Triton X-100	168-11805	化学用	500 mL	9002-93-1
	166-11801		3 L	
	167-21101	生物化学用	3 L	
月桂醇聚氧乙烯醚	162-21075	生物化学用	500 g	9002-92-0
NP-40替代物	142-10061	分子生物学用	100 mL	9016-45-9

稳定剂/保护剂

产品名称	产品编号	等级	规格	CAS号
酪蛋白酸钠	143-02845	化学用	500 g	9005-46-3
聚乙烯吡咯烷酮K30	165-17035	和光特级	500 g	9003-39-8
PVP K30	161-03105	—	500 g	9003-39-8
alpha-环化糊精	033-08332	和光一级	25 g	10016-20-3
beta-环化糊精	034-08345	和光一级	500 g	7585-39-9
阿拉伯树胶	016-00025	—	500 g	9000-01-5
D(-)-甘露醇	133-00845	试剂特级	500 g	69-65-8
L(+)-抗坏血酸钠	196-01252	和光特级	25 g	134-03-2
苯甲酸钠	198-01295	和光特级	500 g	532-32-1

标准品/质控品

产品名称	产品编号	等级	规格	CAS号
D(-)-山梨糖醇	198-03755	和光一级	500 g	50-70-4
L-胱硫醚标准品	030-23801	Traceable Reference Material	50 mg	56-88-2
肌酐	031-04592	和光特级	25 g	60-27-5
	039-04593		10 g	
氯化锂	127-01165	和光特级	500 g	7447-41-8
碳酸锂	126-01135	和光特级	500 g	554-13-2

缓冲体系

产品名称	产品编号	等级	规格	CAS号
咪唑	095-00015	和光特级	500 g	288-32-4
三羟甲基氨基甲烷(Tris碱)	207-06275	试剂特级	500 g	77-86-1
β-甘油磷酸钠水合物	191-02042	—	25 g	819-83-0
	193-02041		100 g	
甘氨酸	073-00732	试剂特级	25 g	56-40-6
	077-00735		500 g	
磷酸二氢铵	019-03315	和光一级	500 g	7722-76-1
磷酸二氢钾	167-24685	药物测试用	500 g	7778-77-0
磷酸二氢钠二水合物	192-02815	试剂特级	500 g	13472-35-0
磷酸氢二钾	164-04295	试剂特级	500 g	7758-11-4
氯化钾	163-03545	试剂特级	500 g	7447-40-7
柠檬酸	030-05525	和光特级	500 g	77-92-9
十水四硼酸钠/硼砂	194-01415	试剂特级	500 g	1303-96-4
碳酸钾	165-03505	和光一级	500 g	584-08-7
碳酸钠	196-01595	和光一级	500 g	497-19-8
三乙醇胺	142-05615	和光一级	500 mL	102-71-6
2-氨基-2-甲基-1-丙醇 (AMP)	014-02586	和光一级	500 mL	124-68-5

底物

产品名称	产品编号	等级	规格	CAS号
2-酮戊二酸	115-00082	和光特级	25 g	328-50-7
L-天冬氨酸	013-04832	试剂特级	25 g	56-84-8
	015-04831		100 g	
	017-04835		500 g	

显色相关

产品名称	产品编号	等级	规格	CAS号
4-氨基氨替吡啉	011-02275	试剂特级	500 g	83-07-8
DA-67	046-22341	生物化学用	50 mg	115871-18-6
DA-64	043-22351	生物化学用	100 mg	115871-19-7
TOOS	346-08811	—	1 g	82692-93-1
还原型辅酶 I 二钠	042-16233	生物化学用	500 mg	606-68-8
	046-16231		100 mg	
链霉亲和素	192-17864	免疫化学用	25 mg	9013-20-1
硫酸铜	030-04442	试剂一级	25 g	7758-98-7
溴甲酚绿	020-02682	试剂特级	25 g	76-60-8

其他

产品名称	产品编号	等级	规格	CAS号
L-精氨酸	011-04615	和光特级	500 g	74-79-3
谷氨酸钠一水合物	198-02035	和光特级	500 g	6106-04-3
6-氨基己酸	012-09645	和光一级	500 g	60-32-2
尿素	219-00175	试剂特级	500 g	57-13-6
柠檬酸钠	191-01785	试剂特级	500 g	6132-04-3
盐酸胍	077-02435	和光特级	500 g	50-01-1

上述试剂仅供实验研究用,不可用作“医药品”、“食品”、“临床诊断”等。

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit our online catalog to search for other products from FUJIFILM Wako: <https://labchem-wako.fujifilm.com> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM

富士胶片 and 光(广州)贸易有限公司

广州市越秀区先烈中路69号东山广场30楼3002、3003、3011室

北京 Tel: 13611333218 上海 Tel: 021 62884751

广州 Tel: 020 87326381 香港 Tel: 852 27999019

询价: wkgz.info@fujifilm.com

官网: labchem.fujifilm-wako.com.cn

官方微信



目录价查询

